

WATERCHECK 1 & 2 SAMPLING INSTRUCTIONS

PLEASE READ ALL THE INSTRUCTIONS COMPLETELY BEFORE COLLECTING YOUR SAMPLE

One test kit represents one sampling location. If you require testing on multiple locations you will need to purchase multiple test kits.

1. Freeze the blue gel ice pack for at least 8 hours prior to sample collection. Return to cooler just prior to sampling and/or return shipping.
2. Remove any aerators, strainers, mixers or purifiers from the tap.
3. Run the cold water for approximately 3-5 minutes, then adjust flow to a thin stream to avoid splashing or overfilling.
4. Fill the bottle to the neck, cap securely and return to the shipping container.
5. Complete the enclosed lab form, **making sure to indicate the date and time you sampled**, type of water (treated/untreated) and sample identification which identifies where you collected your sample (i.e., kitchen sink, well head, etc.). Retain a copy of the lab form for your records and return the completed form to the shipping container.
6. Mail your sample to the following testing location. Sample should be sent via guaranteed overnight delivery.

Testing Location: **National Testing Laboratories, Ltd.
556 South Mansfield Street
Ypsilanti, MI 48197**

7. PLEASE NOTE: Testing takes 5-7 business days from receipt. Results can be emailed or faxed upon completion, please provide your email or fax number on the form provided if not already preprinted.

If you have any questions regarding this Watercheck 1 & 2 Test or any other services we provide, please contact a Technical Service Representatives at 1-800-458-3330 for assistance

Thank you for choosing National Testing Laboratories, Ltd., for your water testing services!